 DE - Sicherheitshinweise - Handzahnbürsten.....	2
 EN - Safety Instructions - Manual Toothbrush	2
 IT - Istruzioni di Sicurezza -Spazzolino manuale	3
 FR - Consignes de Sécurité - brosse à dents manuelle.....	4
 ES - Instrucciones de Seguridad - cepillo de dientes manual	4
 SE - Säkerhetsanvisningar - manuell tandborste.....	5
 PL - Instrukcje Bezpieczeństwa - manualna szczoteczka do zębów	6
 NL - Veiligheidsinstructies - handtandenborstel	7

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EN - Safety Instructions - Manual Toothbrush

General Safety Guidelines for Using a Manual Toothbrush

- ▶ **Avoid Excessive Pressure:** Brushing your teeth with too much pressure can damage your gums and wear down enamel. Brush gently using circular motions.
- ▶ **Replace Your Toothbrush Regularly:** Toothbrushes should be replaced every 2–3 months or sooner if the bristles are worn out. Worn bristles are less effective at cleaning teeth.
- ▶ **Rinse Thoroughly After Use:** After brushing, rinse the toothbrush thoroughly with water and store it in a dry place to prevent bacterial growth.
- ▶ **Do Not Share Your Toothbrush:** Sharing toothbrushes can spread germs or infections and should be avoided.

- ▶ **Use Soft Bristles:** Opt for a toothbrush with soft bristles to protect your gums and enamel. Hard bristles can cause gum irritation and enamel damage.
- ▶ **Apply Proper Brushing Technique:** Use circular motions while brushing and ensure all tooth surfaces are cleaned for thorough oral hygiene.
- ▶ **Store Upright:** Avoid placing the toothbrush on surfaces that may be contaminated with germs after use.
- ▶ **Avoid Surface Contact:** If you have open wounds or injuries in your gums, avoid using the toothbrush until the area has healed.
- ▶ **Monitor Gum Bleeding:** Light gum bleeding may be normal at first, but if it persists, worsens, or your gums feel sensitive, consult a dentist.
- ▶ **Choose Toothbrushes According to Age and Needs:** Select toothbrushes that meet specific needs (e.g., special models for children, sensitive teeth, or braces) for optimal results and protection.

⚠ Always refer to the product packaging or instructions for specific recommendations.

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