 DE - Sicherheitshinweise - Zahnseide	2
 EN - Safety Instructions - Dental floss	2
 IT - Istruzioni di Sicurezza -Filo interdentale.....	3
 FR - Consignes de Sécurité - Fil dentaire.....	4
 ES - Instrucciones de Seguridad - Hilo dental.....	5
 SE - Säkerhetsanvisningar - Tandtråd.....	6
 PL - Instrukcje Bezpieczeństwa - Nić dentystyczna	6
 NL - Veiligheidsinstructies - Tandfloss	7

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EN - Safety Instructions - Dental floss

⚠ General Safety Warnings for Using Dental Floss ⚠

- ▶ **Gentle application:** Dental floss should be gently guided between the teeth to avoid gum injuries. Excessive pressure can cause gum bleeding and irritation.

- ▶ Do not insert too deeply into the gums: Ensure the floss does not go too far under the gums, as this can damage the tissue.
- ▶ Use only once: Each piece of floss should be used only once to prevent the spread of bacteria.
- ▶ Regular use: Use dental floss daily to remove plaque between teeth and reduce the risk of gum disease and cavities. Light gum bleeding may be normal at the start of flossing, but if bleeding persists or worsens, consult a dentist.
- ▶ Store floss safely: Floss is a supplement to daily oral care but does not replace thorough tooth brushing.
- ▶ Not a substitute for brushing: Dental floss complements oral hygiene but does not replace tooth brushing or mouthwash. It should only be used as an additional step in daily oral care.
- ▶ Caution with bridges and implants: Special floss or interdental brushes should be used for bridges, implants, or braces to avoid damage.
- ▶ Floss for sensitive gums: People with sensitive or inflamed gums should use soft or waxed floss to minimize irritation.
- ▶ Use proper technique: To use floss effectively, curve it gently into a C-shape around the tooth and move it carefully up and down to remove plaque.
- ▶ Regular use: Daily flossing is essential to remove plaque between teeth and lower the risk of gum disease and cavities.
- ▶ For children aged 6 and older: Once children have tightly spaced teeth, around age 6, they should begin flossing (using floss picks or floss holders), as toothbrushes often cannot reach all spaces. Children should initially floss under the guidance and supervision of an adult.

 **Always follow the instructions on the product or packaging.**

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